

SYRINGE SHIELDS

ALL CONFIGURATIONS



**NUCLEAR
SHIELDS**



NUCLEAR SHIELDS



SYRINGE SHIELDS

Our tungsten syringe shields designed to ensure optimized attenuation and reduce hand and finger doses. These syringe shields are shielded with 2,1 mm tungsten and a 5.2 density lead glass viewing window. The body of the tungsten syringe shields is finished with anodized aluminum for easy cleaning and decontamination. These tungsten syringe shields are designed for 1-10 ml syringes. Most syringe shield volumes are available in multiple internal dimensions for different syringe manufacturers.

Big Screw Product Code	Clip-Lock Product Code	Volume	inside dimensions
SS2BS1-1	SS2CL1-1	1 ml	8.8 x 72 mm
SS2BS1-2	SS2CL1-2	1 ml	8.8 x 76 mm
SS2BS1-3	SS2CL1-3	1 ml	11.5 x 76 mm
SS2BS2-4	SS2CL2-4	2 ml	11.5 x 53 mm
SS2BS2-5	SS2CL2-5	2 ml	11.5 x 60 mm
SS2BS3-6	SS2CL3-6	3 ml	11.5 x 65 mm
SS2BS3-7	SS2CL3-7	3 ml	11.5 x 69 mm
SS2BS5-8	SS2CL5-8	5 ml	15.4 x 61 mm
SS2BS5-9	SS2CL5-9	5 ml	15.4 x 65 mm
SS2BS5-10	SS2CL5-10	5 ml	15.4 x 68 mm
SS2BS10-11	SS2CL10-11	10 ml	18.2 x 75 mm
SS2BS10-12	SS2CL10-12	10 ml	18.2 x 85 mm



NUCLEAR SHIELDS



LEAD GLASS SYRINGE SHIELDS

Our syringe shields with lead glass viewing window for optimal visibility of the syringe. The clip-lock mechanism will make sure that the syringe is held in place. The anodized aluminum end cap is lined with 2 mm tungsten shielding on the inside to also provide shielding at the clip-lock mechanism, where the lead glass is not able to reach. This makes sure that you are optimally protected. The lead glass has been designed with rounded edges to increase the strength.

Lead Glass Product Code	Volume	Inside dimensions
SS2LG1-1	1 ml	8.8 x 72 mm
SS2LG1-2	1 ml	8.8 x 76 mm
SS2LG1-3	1 ml	11.5 x 76 mm
SS2LG2-4	2 ml	11.5 x 53 mm
SS2LG2-5	2 ml	11.5 x 60 mm
SS2LG3-6	3 ml	11.5 x 65 mm
SS2LG3-7	3 ml	11.5 x 69 mm
SS2LG5-8	5 ml	15.4 x 61 mm
SS2LG5-9	5 ml	15.4 x 65 mm
SS2LG5-10	5 ml	15.4 x 68 mm
SS2LG10-11	10 ml	18.2 x 75 mm
SS2LG10-12	10 ml	18.2 x 85 mm



NUCLEAR SHIELDS



PET SYRINGE SHIELDS

Our tungsten syringe shields for PET with 8mm tungsten and our convenient Big-Screw locking mechanism. They are available in multiple sizes to match your required barrel dimensions. Please refer to the size chart below, or download the PET syringe shield size chart document below.

PET Product Code	Volume	Inside dimensions
SS2BS1-1P	1 ml	8.8 x 72 mm
SS2BS1-2P	1 ml	8.8 x 76 mm
SS2BS1-3P	1 ml	11.5 x 76 mm
SS2BS2-4P	2 ml	11.5 x 53 mm
SS2BS2-5P	2 ml	11.5 x 60 mm
SS2BS3-6P	3 ml	11.5 x 65 mm
SS2BS3-7P	3 ml	11.5 x 69 mm
SS2BS5-8P	5 ml	15.4 x 61 mm
SS2BS5-9P	5 ml	15.4 x 65 mm
SS2BS5-10P	5 ml	15.4 x 68 mm
SS2BS10-11P	10 ml	18.2 x 75 mm
SS2BS10-12P	10 ml	18.2 x 85 mm



CONTACT FORM